

# Sterile Processing Technician Practice Exam (Sample)

## Study Guide



**Everything you need from our exam experts!**

**Copyright © 2026 by Examzify - A Kaluba Technologies Inc. product.**

**ALL RIGHTS RESERVED.**

**No part of this book may be reproduced or transferred in any form or by any means, graphic, electronic, or mechanical, including photocopying, recording, web distribution, taping, or by any information storage retrieval system, without the written permission of the author.**

**Notice: Examzify makes every reasonable effort to obtain accurate, complete, and timely information about this product from reliable sources.**

**SAMPLE**

# Table of Contents

**Copyright** ..... 1

**Table of Contents** ..... 2

**Introduction** ..... 3

**How to Use This Guide** ..... 4

**Questions** ..... 5

**Answers** ..... 8

**Explanations** ..... 10

**Next Steps** ..... 16

SAMPLE

# Introduction

Preparing for a certification exam can feel overwhelming, but with the right tools, it becomes an opportunity to build confidence, sharpen your skills, and move one step closer to your goals. At Examzify, we believe that effective exam preparation isn't just about memorization, it's about understanding the material, identifying knowledge gaps, and building the test-taking strategies that lead to success.

This guide was designed to help you do exactly that.

Whether you're preparing for a licensing exam, professional certification, or entry-level qualification, this book offers structured practice to reinforce key concepts. You'll find a wide range of multiple-choice questions, each followed by clear explanations to help you understand not just the right answer, but why it's correct.

The content in this guide is based on real-world exam objectives and aligned with the types of questions and topics commonly found on official tests. It's ideal for learners who want to:

- Practice answering questions under realistic conditions,
- Improve accuracy and speed,
- Review explanations to strengthen weak areas, and
- Approach the exam with greater confidence.

We recommend using this book not as a stand-alone study tool, but alongside other resources like flashcards, textbooks, or hands-on training. For best results, we recommend working through each question, reflecting on the explanation provided, and revisiting the topics that challenge you most.

**Remember:** successful test preparation isn't about getting every question right the first time, it's about learning from your mistakes and improving over time. Stay focused, trust the process, and know that every page you turn brings you closer to success.

Let's begin.

# How to Use This Guide

**This guide is designed to help you study more effectively and approach your exam with confidence. Whether you're reviewing for the first time or doing a final refresh, here's how to get the most out of your Examzify study guide:**

## **1. Start with a Diagnostic Review**

**Skim through the questions to get a sense of what you know and what you need to focus on. Your goal is to identify knowledge gaps early.**

## **2. Study in Short, Focused Sessions**

**Break your study time into manageable blocks (e.g. 30 - 45 minutes). Review a handful of questions, reflect on the explanations.**

## **3. Learn from the Explanations**

**After answering a question, always read the explanation, even if you got it right. It reinforces key points, corrects misunderstandings, and teaches subtle distinctions between similar answers.**

## **4. Track Your Progress**

**Use bookmarks or notes (if reading digitally) to mark difficult questions. Revisit these regularly and track improvements over time.**

## **5. Simulate the Real Exam**

**Once you're comfortable, try taking a full set of questions without pausing. Set a timer and simulate test-day conditions to build confidence and time management skills.**

## **6. Repeat and Review**

**Don't just study once, repetition builds retention. Re-attempt questions after a few days and revisit explanations to reinforce learning. Pair this guide with other Examzify tools like flashcards, and digital practice tests to strengthen your preparation across formats.**

**There's no single right way to study, but consistent, thoughtful effort always wins. Use this guide flexibly, adapt the tips above to fit your pace and learning style. You've got this!**

## Questions

SAMPLE

- 1. All things can cause cycle cancellation in LTGP Sterilization except?**
  - A. Tyvek**
  - B. Paper - plastic pouch**
  - C. Cellulose based materials**
  - D. Moisture**
  
- 2. What should be done with equipment that has a cracked or broken electrical cord?**
  - A. Cleaned then repaired**
  - B. Discarded**
  - C. Sent out for repair**
  - D. Cleaned, disinfected, and sent out for repair**
  
- 3. What is the only organ system that can be surgically removed without impeding one's activities of daily living?**
  - A. Reproductive System**
  - B. Respiratory System**
  - C. Urinary System**
  - D. Integumentary System**
  
- 4. What is the coldest part of a steam sterilizer?**
  - A. Drain line**
  - B. Jacket**
  - C. Chamber**
  - D. Gasket**
  
- 5. Which parameter is not directly measured in a steam sterilization process?**
  - A. Temperature**
  - B. Pressure**
  - C. Humidity**
  - D. Time**

- 6. Which of the following is a critical step in the reprocessing of medical instruments?**
- A. Cleaning**
  - B. Disinfecting**
  - C. Sterilization**
  - D. All of the above**
- 7. Which body poses the standard for surgical instruments?**
- A. FDA**
  - B. CDC**
  - C. OSHA**
  - D. AAMI**
- 8. What type of fire extinguisher is appropriate for use on combustible metals?**
- A. Type K**
  - B. Type E**
  - C. Type B**
  - D. Type D**
- 9. Which of the following circumstances can prevent product degradation?**
- A. Avoiding temperature and humidity extremes**
  - B. Storing instruments directly on the floor**
  - C. Storing instruments within 20" from the sprinkler system**
  - D. Both a and c**
- 10. Why is it important to keep sterilized items in a controlled environment?**
- A. To ensure easy access for staff**
  - B. To maintain cleanliness and prevent contamination**
  - C. To allow for increased airflow**
  - D. To reduce packaging costs**

## Answers

SAMPLE

1. A
2. D
3. A
4. A
5. C
6. D
7. D
8. D
9. D
10. B

SAMPLE

## **Explanations**

SAMPLE

**1. All things can cause cycle cancellation in LTGP Sterilization except?**

- A. Tyvek**
- B. Paper - plastic pouch**
- C. Cellulose based materials**
- D. Moisture**

In the context of Low Temperature Gas Plasma (LTGP) sterilization, understanding the materials that can impact the sterilization cycle is crucial for the effective use of this method. Tyvek is often recognized as a compatible material for LTGP sterilization because it allows the gas to penetrate and effectively sterilize the instruments or devices it encases. This permeability is essential for ensuring that the sterilization agent can contact the items inside, which helps in achieving proper sterilization. In contrast, other materials mentioned—such as paper-plastic pouches, cellulose-based materials, and moisture—can interfere with the LTGP process. Paper-plastic pouches may hinder gas penetration due to their material properties, cellulose-based materials may absorb moisture or impede gas movement, and excessive moisture is generally known to cause cycle cancellations since it can lead to ineffective sterilization. Thus, Tyvek's unique properties make it the only item listed that does not contribute to the cancellation of the sterilization cycle, allowing it to be safely used in LTGP processes without compromising effectiveness.

**2. What should be done with equipment that has a cracked or broken electrical cord?**

- A. Cleaned then repaired**
- B. Discarded**
- C. Sent out for repair**
- D. Cleaned, disinfected, and sent out for repair**

When dealing with equipment that has a cracked or broken electrical cord, the most appropriate course of action is to clean, disinfect, and send it out for repair. This approach prioritizes both safety and equipment functionality. Firstly, cleaning the equipment ensures that any contaminants on its surface are removed. This step is vital in maintaining a safe working environment, especially in settings like sterile processing, where hygiene is critical. Following the cleaning process, disinfecting is necessary to eliminate any remaining microorganisms, further ensuring the safety of personnel and patients. Sending the damaged equipment out for repair is essential because it allows for a qualified technician to address the electrical fault properly. Using equipment with a compromised electrical cord poses a significant safety risk, including potential electrical shocks or equipment failure. In this scenario, simply discarding the equipment without attempting to repair or maintain it would be wasteful, particularly if the equipment is otherwise functional. Repairing rather than discarding also aligns with sustainability practices and responsible resource management within healthcare settings.

**3. What is the only organ system that can be surgically removed without impeding one's activities of daily living?**

- A. Reproductive System**
- B. Respiratory System**
- C. Urinary System**
- D. Integumentary System**

The correct answer is the reproductive system because it is the only organ system that can be completely surgically removed without fundamentally affecting a person's ability to carry out daily activities. While the removal of reproductive organs can have significant emotional and psychological impacts, the primary functions necessary for survival and daily functioning (e.g., digestion, respiration, circulation) are not directly dependent on the reproductive system. For instance, individuals can still perform essential life activities, such as eating, walking, and interacting with others, even after the surgical removal of reproductive organs, be it in a male or female context. Hence, while there may be hormonal changes and potential impacts on quality of life, daily living activities are largely unaffected. In contrast, the other organ systems have crucial roles that directly relate to survival and daily functioning. The respiratory system is essential for gas exchange and breathing, the urinary system plays a vital role in waste elimination, and the integumentary system (which includes skin, hair, and nails) provides barriers, protection, and sensory reception. The removal of any of these systems would severely impede daily living and require significant adjustments.

**4. What is the coldest part of a steam sterilizer?**

- A. Drain line**
- B. Jacket**
- C. Chamber**
- D. Gasket**

The coldest part of a steam sterilizer is typically found within the drain line. In a steam sterilization process, the drain line is where condensate—a byproduct of steam—accumulates and exits the sterilizer. Since this line is not exposed to the continuous flow of high-temperature steam like other areas of the sterilizer, it tends to cool down more quickly, making it cooler compared to other parts of the system. Understanding this is crucial for sterilization processes, as it highlights the importance of ensuring that all items being processed reach the required temperature for effective sterilization. Proper functioning of the drain line is also essential to maintain the overall efficiency of the sterilizer and prevent any potential cooling of steam that could impact the sterilization cycle. The jacket, chamber, and gasket each have their roles, but they are more directly involved with the heating process or maintaining the seal required for sterilization. The jacket surrounds the chamber to create additional heat insulation, while the chamber is the area where the items are sterilized, and the gasket serves to maintain an airtight seal during the cycle.

**5. Which parameter is not directly measured in a steam sterilization process?**

- A. Temperature**
- B. Pressure**
- C. Humidity**
- D. Time**

In a steam sterilization process, the parameters that are directly measured include temperature, pressure, and time, all of which are critical to ensuring effective sterilization. Temperature is essential because sufficient heat is required to kill microorganisms, while pressure is necessary to maintain the temperature of the steam and prevent its condensation. Time is also a crucial factor since prolonged exposure to the sterilizing conditions ensures that all microorganisms are effectively eradicated. Humidity, while it can play a role in the effectiveness of sterilization by influencing steam quality, is not a parameter that is directly measured during the process itself. Instead, it is often inferred from conditions such as temperature and pressure. The focus in sterilization monitoring is primarily on the direct parameters that directly indicate that effective sterilization has occurred, making humidity less central to the immediate measurement strategy during steam sterilization.

**6. Which of the following is a critical step in the reprocessing of medical instruments?**

- A. Cleaning**
- B. Disinfecting**
- C. Sterilization**
- D. All of the above**

In the context of reprocessing medical instruments, each step contributes significantly to ensuring that instruments are safe for patient use, but all these steps work collectively to achieve comprehensive reprocessing. Cleaning is the foundational step, where visible soil, blood, and other contaminants are removed from instruments. If this step is skipped or inadequately performed, it compromises the effectiveness of subsequent processes like disinfection and sterilization, as these processes cannot work properly on dirty instruments. Disinfecting follows cleaning, targeting pathogens that may remain after the initial cleaning. This step reduces the number of viable microorganisms, preparing instruments for the final step. Sterilization is the ultimate goal in reprocessing, where all forms of microbial life, including spores, are killed. This step is vital for items that come into contact with sterile body areas or enter sterile tissues and is critical for patient safety. By understanding that all these processes—cleaning, disinfecting, and sterilization—are essential and interdependent, it becomes evident that they collectively represent the critical steps in the reprocessing of medical instruments. Therefore, recognizing that all these steps are necessary illustrates the complex nature of infection control in healthcare settings.

## 7. Which body poses the standard for surgical instruments?

- A. FDA
- B. CDC
- C. OSHA
- D. AAMI**

The Association for the Advancement of Medical Instrumentation (AAMI) plays a crucial role in setting the standards for surgical instruments. AAMI works to develop and publish guidelines that ensure the safety and effectiveness of medical instruments, particularly those used in healthcare settings. Their standards focus on aspects such as reprocessing, sterilization methods, and maintenance of surgical instruments, which directly impact patient safety and outcomes. AAMI's standards are developed through a consensus of professionals in the field, incorporating the latest research and best practices. By adhering to AAMI standards, sterile processing technicians can guarantee that instruments are properly cleaned, sterilized, and maintained, thereby reducing the risk of infection and ensuring high-quality patient care. In contrast, while the FDA regulates medical devices and oversees their safety and efficacy, it does not set specific standards for surgical instrument reprocessing. The CDC provides guidelines for infection control in healthcare settings, but it does not inherently specify standards for the instruments themselves. OSHA focuses on occupational safety and health, covering workplace safety regulations rather than technical standards specific to surgical instruments. Thus, AAMI stands out as the authority specifically geared toward the standards of surgical instruments.

## 8. What type of fire extinguisher is appropriate for use on combustible metals?

- A. Type K
- B. Type E
- C. Type B
- D. Type D**

The appropriate type of fire extinguisher for use on combustible metals is the Type D extinguisher. This classification is specifically designed to combat fires involving metals such as magnesium, titanium, and sodium, which require unique extinguishing agents due to their reactive properties when exposed to water or common extinguishing agents. Type D extinguishers often contain dry powder agents that work by smothering the fire and interrupting the combustion process. These agents are capable of handling the intense and specific types of fires that combustible metals can produce, which is essential in environments such as manufacturing or laboratories where these materials are present. In contrast, other fire extinguisher types target different fire classes. For example, a Type K extinguisher is used for fires involving cooking oils and fats, Type E is not commonly recognized in the standard classifications (often it might refer to electrical fires in some contexts, but it is not universally accepted), and Type B extinguishers are meant for flammable liquids and gases, which wouldn't be effective on combustible metals. Therefore, understanding the specific characteristics of fire extinguishers based on the materials involved in a fire is vital for safe and effective fire management.

**9. Which of the following circumstances can prevent product degradation?**

- A. Avoiding temperature and humidity extremes**
- B. Storing instruments directly on the floor**
- C. Storing instruments within 20" from the sprinkler system**
- D. Both a and c**

Preventing product degradation is critical in maintaining the integrity and effectiveness of sterile instruments and equipment. Avoiding temperature and humidity extremes directly contributes to preserving the quality and functionality of these items. High temperatures can lead to chemical breakdown, while excessive humidity can promote the growth of microorganisms and cause corrosion. Therefore, maintaining a stable environment minimizes the risk of deterioration. Storing instruments within 20 inches from the sprinkler system, although it may seem counterintuitive, is also relevant. Sprinkler systems are designed to activate during emergencies, and while instruments should not be directly exposed to potential water damage, being within the necessary range ensures compliance with safety codes, which can help in maintaining overall preservation protocols in a sterile processing area. By avoiding extremes in temperature and humidity while adhering to proper storage guidelines relative to fire safety equipment, the chances of product degradation are significantly reduced, ensuring that instruments remain sterile and functional when needed.

**10. Why is it important to keep sterilized items in a controlled environment?**

- A. To ensure easy access for staff**
- B. To maintain cleanliness and prevent contamination**
- C. To allow for increased airflow**
- D. To reduce packaging costs**

Maintaining sterilized items in a controlled environment is essential primarily to uphold their sterility. The main goal of sterilizing instruments and supplies is to eliminate all forms of microbial life, including bacteria, viruses, and spores. Once items are sterilized, exposure to air, moisture, or unfiltered handling can risk contamination. A controlled environment helps regulate factors such as temperature, humidity, and cleanliness, which are crucial in preserving the integrity of the sterility. A sterile item stored in a clean, dry environment is less likely to come into contact with pollutants or pathogens. Additionally, storage conditions can affect the efficacy of the packaging used for sterilized items; if the packaging is compromised, it may no longer provide the necessary protection against contaminants. By ensuring that sterilized items remain in a controlled environment, facilities can better guarantee that these items will remain safe for patient use when needed, thus reducing the risk of hospital-acquired infections or complications associated with surgical procedures.

## Next Steps

**Congratulations on reaching the final section of this guide. You've taken a meaningful step toward passing your certification exam and advancing your career.**

**As you continue preparing, remember that consistent practice, review, and self-reflection are key to success. Make time to revisit difficult topics, simulate exam conditions, and track your progress along the way.**

**If you need help, have suggestions, or want to share feedback, we'd love to hear from you. Reach out to our team at [hello@examzify.com](mailto:hello@examzify.com).**

**Or visit your dedicated course page for more study tools and resources:**

**<https://sterileprocessingtech.examzify.com>**

**We wish you the very best on your exam journey. You've got this!**

SAMPLE